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owner's manual

RAIOS X

X-RAY

Spectro 70X

*Seletronic
Eletronic*

DABI ATLANTE

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Raio X Spectro 70X


fig. A


fig. B


fig. C


fig. D

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Contents

Legend	18
Installation	19
Procedures prior to equipment use/ reutilization	19
Language selection (Seletronic)	19
Operation	19
Exposure time adjustment (Seletronic)	19
Standard time table	20
Radiographic techniques	20
Programming (Seletronic) 16	
Standard data restoration (Seletronic)	20
Cleaning	20
Maintenance	20
Recommendations and care	20
Special handling precautions for the movable column version	21
Points to be observed prior to the first load application	21
Method used for measuring the X-ray tube voltage	21
Method used for measuring the X-ray tube current	21
Procedure in case equipment functioning is altered	21
Sensitivity to predictable environmental conditions in normal working situations	21
Procedures in case of equipment disposal	21
Transportation and storage conditions	21
After installation	22
Preparation prior to delivery	22
Dabi Atlante authorized service network ...	22
Identification labels	22
Symbols	23
Final recommendations	23
Equipment warranty	23
Technical features	24
Troubleshooting	25
Appendices	26

Presentation

Dear Customer,

Welcome!

We are really proud of your preference for our products. Your familiarization with the **Spectro X Dental X-Ray Machine** will help maximize its performance and increase its durability and benefits. Therefore, before starting your equipment, be sure to read this manual very attentively.

Important

All the information, illustrations and specifications contained in this manual are based on data existing at the time of its publication. We reserve the right to make any modifications deemed necessary both on the product and in this manual without previous notice.

Equipment description

This equipment has been designed for controlled emission of ionizing radiation with the purpose to produce radiographic images for dental diagnostic and treatment procedures using appropriate film types. It features all the required resources to assist dental professionals in the maintenance of high-quality imaging so as to enable proper assessment and diagnosis of the area to be treated.

Except for the "Seletronic" and "Eletronic" sections, all other information contained in this manual is applicable to the **Spectro 70X Dental X-Ray Machine models**.

Legend

Spectro 70X Seletronic (fig. A)

- ① X-ray tubehead
- ② Arms
- ③ Timer box
- ④ Equipment protection fuse
- ⑤ Master switch
- ⑥ Remote control (Seletronic or Eletronic)
- ⑦ Remote control cable (spiraled)
- ⑧ Post
- ⑨ Supporting base
- ⑩ Caster
- ⑪ Focal point: indicates the exact position of the x-ray emitting center.
- ⑫ Angulation indicator: graduated from 10 to 10 degrees, indicates the head's inclination angle
- ⑬ Primary collimator: it limits the x-ray beam to 5.9 cm in diameter at the locating cylinder outlet
- ⑭ Locating cylinder: it positions the film within the x-ray emitting area, ensuring the focus/skin distance will not be shorter than 20 cm shorter. It has a secondary collimator to eliminate shades in the radiograph
- (A) (B) (C) (D) Identification tags
- ⑮ ⑯ ⑰ ⑱ Optional items only for the Eletronic model

Control legend (Seletronic) (fig. B)

- ① Display
- ② Exposure time increase key
- ③ Exposure time decrease key
- ④ Ionizing radiation emission indicator
- ⑤ Triggering key
- ⑥ Save key
- ⑦ Restore key
- ⑧ Anterior interproximal key
- ⑨ Posterior interproximal key
- ⑩ Periapical keys
- ⑪ Jaw-occlusal key
- ⑫ Maxilla-occlusal key
- ⑬ Film type selection key
- ⑭ Adult/child selection key

Display Indicators (Seletronic) (fig. C)

- ① Indicates the type of radiograph available:
PER= Periapical
INTERPROX. ANT. = Anterior Interproximal
INTERPROX. POS. = Posterior interproximal

OCCLUSAL MAXILAR = Maxillary occlusal

OCCLUSAL MANDIB. = Mandibular occlusal

- ② Indicates the type of film available for the function: f:D or f:E or f:F or F:DIG (digital)
- ③ Indicates that either the child or adult option is available for the function. Ad = Adult or INF = Child
- ④ Indicates the exposure time available for the selected function.
- ⑤ Indicates that the standard time was modified by the user.

Spectro 70X Eletronic (fig. A)

- ① X-ray tubehead
- ② Arms
- ③ Timer box
- ④ Equipment protection fuse
- ⑤ Master switch
- ⑥ Remote control (Seletronic or Eletronic)
- ⑦ Remote control cable (spiraled)
- ⑧ Post
- ⑨ Supporting base
- ⑩ Caster
- ⑪ Focal point: indicates the exact position of the x-ray emitting center.
- ⑫ Angulation indicator: graduated from 10 to 10 degrees, indicates the head's inclination angle.
- ⑬ Primary collimator: it limits the x-ray beam to 5.9 cm in diameter at the locating cylinder outlet.
- ⑭ Locating cylinder: it positions the film within the x-ray emitting area, ensuring the focus/skin distance will not be shorter than 20 cm shorter. It has a secondary collimator to eliminate shades in the radiograph.
- ⑮ Voltage adjustment selection switch (optional item - Eletronic)
- ⑯ High-voltage indication LED (RED) (optional item - Eletronic)
- ⑰ Nominal-voltage indication LED (GREEN) (optional item - Eletronic)
- ⑱ Low-voltage indication LED (RED) (optional item - Eletronic)
- (A) (B) (C) (D) Identification tags.

Control (Eletronic) (fig. B)

- ① Display
- ② Exposure time increase key
- ③ Exposure time decrease key
- ④ Ionizing radiation emission indicator
- ⑤ Triggering key

Under and Overload Indicators (Eletronic) (fig. D)

EUL – Indicates that the mains supply voltage is low and that the equipment will not shoot.

EUH – Indicates that the mains supply voltage is high and that the equipment will not shoot.

X-radiation Source Set (fig. F)

- ① Head
- ② X-ray tube
- ③ Oil
- ④ Acrylic
- ⑤ Aluminum filter (see table "quality equivalency filtration")
- ⑥ Primary collimator (diaphragm)
- ⑦ Secondary Collimator (diaphragm)
- ⑧ Sealing cover
- ⑨ Collimator cylinder
- ⑩ Focal point

Installation

This equipment must be unpacked and installed by an authorized Dabi Atlante technician; otherwise, warranty will be voided.

Procedures prior to equipment use/reutilization

Prior to equipment use/reutilization, follow the operation and cleaning instructions in this manual.

Language selection (Seletronic)

When starting the machine, keep the "Adult/Child" key pressed.

When the screen opens, the display will show three languages (Portuguese, English and Spanish), which can be changed by the "+" and "-" keys.

Press the "S" key to save the language chosen.

Operation**Seletronic**

Position the locating cylinder according to the radiographic taking type desired.

Turn on the master switch. The message "DABIATLANTE VERSION 1.1" will show on the display, which will be followed by the message "SELECT RADIOGRAPH".

Select the radiograph type using the anterior Interproximal, posterior Interproximal, Occlusal-Lower Jaw, Occlusal-Upper Jaw and Periapical keys.

Note that the periapical function is selected by using the control keys which reproduce the dental arch.

Press the "adult/child" key to select the adult or child option.

Select the film type to be used.

Stay away at a safe distance (see item: Recommendations, care and warnings).

Keep the shooting key pressed to activate the equipment. After a sound warning is heard the display will show "REAL TIME __ _ s COMPENS. +/- __ _ S".

After shooting, the display will show the actual exposure time as the Timer features an intelligent system that compensates the x-ray application time according to variations in the mains voltage. For this reason, the actual time described on the display after shooting must be regarded as the exposure time.

After use, turn off the equipment master switch.

Releasing the shooting key before the selected time has elapsed interrupts x-ray emission immediately and the message "REAL TIME __ _ s SHOOT INTERRUPT." will show on the display, indicating the actual shooting time.

Teeth 16/17/18, 26/27/28, 36/37/38, 46/47/48 cannot be selected in the child option.

The time selected by the operator for periapical radiographs must be similar for both arches.

Eletronic

Position the locating cylinder according to the radiographic taking type desired.

Turn on the master switch.

In case the equipment has a voltage regulator (optional item), turn the voltage selector to adjust voltage. Turn it until the central green LED in the timer case goes on. While a red LED remains on, the equipment will be disabled.

Select the exposure time by plus key (+) and minus key (-) to the following values : 0.2; 0.25; 0.32; 0.4; 0.5; 0.63; 0.8; 1.0; 1.25; 1.6; 2.0; 2.5; 3.2 seconds (see item "Radiographic Techniques" in this manual).

Stay away at a safe distance (see item: Recommendations, care and warnings).

Keep the shooting key pressed to activate the equipment. The yellow LED will go on and sound warning will be heard, indicating that the equipment will be emitting radiation during the selected time.

Releasing the shutter button before the end of the selected time, stop the emission of X-ray immediately and the display will show the "IE" message (Interrupted Error).

Wait a few seconds before shooting again. This is the cooling time required which is indicated by the blinking signal on the display. During such time, the keys will be inactive and will show a countdown of the cooling time on the display.

After use, turn off the equipment master switch.

To obtain the load application time, 0,360 must be added (ampoule pre-warming time) to any value on the display.

Exposure time adjustment (Seletronic)

After informing: radiograph type; adult/child and film type, the Seletronic automatically determines

the exposure time, which can be extended/reduced by using keys “+” or “-”.

In case the operator changes the standard time, an asterisk will show next to the selected time and will remain next to that value when it is saved.

Standard time table

Periapical-Teeth	E Film	D Film	F Film	Digital
11 to 13 21 to 23	0,50	1,00	0,30	0,25
14-15 24-25	0,55	1,10	0,35	0,30
16-17-18 26 27 28	0,60	1,20	0,35	0,30
31 to 33 41 to 43	0,45	0,90	0,25	0,20
34 to 38 44 to 48	0,50	1,00	0,30	0,25
51 to 53 61 63	0,35	0,70	0,20	0,15
54 55 64 65	0,40	0,80	0,25	0,20
71 to 73 81 to 83	0,30	0,60	0,20	0,15
84 85 74 75	0,35	0,70	0,20	0,15

Adult

Radiograph Type	E Film	D Film	F Film	Digital
Interprox. Anterior	0,50	1,00	0,30	0,25
Interprox. Posterior	0,60	1,20	0,35	0,30
Occlusal Maxillary	0,80	1,60	0,50	0,40
Occlusal Mandibular	0,70	1,40	0,40	0,35

Infantil

Radiograph Type	E Film	D Film	F Film	Digital
Interprox. Anterior	0,35	0,70	0,20	0,15
Interprox. Posterior	0,40	0,80	0,25	0,20
Occlusal Maxillary	0,60	1,20	0,35	0,30
Occlusal Mandibular	0,50	1,00	0,30	0,25

Radiographic techniques

Radiographic techniques that use a focus/film distance of 25 cm

Technique	Region	Group D	Group E
Bisectrix periapical (“short cone”)	Maxilla	0,4 to 0,5	0,3 to 0,4
	Jaw	0,3 to 0,4	0,2 to 0,3
Interproximal (“bite wing”)	Posterior teeth	0,4 to 0,5	0,3 to 0,4
Parallelism periapical (“long cone”)	Maxilla	0,6 to 0,8	0,5 to 0,6
	Jaw	0,5 to 0,6	0,4 to 0,5

Radiographic techniques that use a focus/film distance of 30 cm

Technique	Region	Group D	Group E
Occlusal	Maxilla	0,6 to 0,8	0,5 to 0,6
	jaw	0,5 to 0,6	0,4 to 0,5

Programming (Seletronic)

Select and set radiograph type (see item “Operation” in this Manual).

Set exposure time (see item “Exposure Time Adjustment” in this Manual).

Press key “S”. The message “SAVE CHANGED DATA?” will show on the display.

Press key “S” again to save the setting or press any other key to cancel the saving process. The message “SAVING FILM + TIME” will show on the display.

Standard data restoration (Seletronic)

Press key “R” to restore the standard time for the selected function. The message “RESTORE STANDARD DATA?” will show on the display.

Press key “R” again to restore the standard setting or press any other key to cancel. The message “RESTORING STANDARD TIME” will show on the display.

In case the operator presses any other key, the message “CANCELLED” will show on the display.

Cleaning

Clean the equipment surface with a clean soft cloth moistened in Aplic Odonto by Dabi Atlante or another product with similar chemical properties.

The use of other chemical products is not recommended as it may damage the equipment.

Maintenance

Preventive maintenance

In order to reduce equipment faults and increase durability, call Dabi Atlante Authorized Technical Assistance and set a regular preventive maintenance plan.

Corrective maintenance

In case the equipment presents any faults that are not listed in “Troubleshooting”, contact Dabi Atlante Authorized Technical Assistance.

Do not open the equipment and/or try to repair it as such action may aggravate the problem or even generate further malfunctioning.

Circuit diagrams, list of parts, components, others

In case they are necessary, ask Dabi Atlante for electric diagrams and/or lists of parts, components or other information, which will be supplied by agreement.

Recommendations and care

The requirements must be met and the equipment must be installed and operated

according to EMC information provided by the attached documents.

Portable and mobile radiofrequency communication equipment can affect the Dental X-Ray Spectro 70X models Seletronic and Eletronic

This equipment must be exclusively used by health care professionals as it may cause radio interference or interrupt the operation of other nearby equipment. Mitigatory easures, such as re-orientation, re-location of equipment or shielding of the site, may be necessary.

The operator must wear a lead coat and vest during exposure. A lead coat or vest are not supplied with the equipment.

During exposure, the operator must keep a minimum distance of 2m from the patient's head or stand behind a physical barrier.

Install the equipment in an appropriate site which is protected from sunlight and humidity.

Unplug the equipment or turn off the clinic's general switch when it is not used for long periods of time.

Check the equipment's voltage before plugging it in.

This equipment is meant for dental use only and must be operated by qualified personnel according to instructions in this Manual.

In case the recommendations and precautions concerning this equipment are not observed, side effects may affect its operators and patients due to the presence of ionizing radiation.

Do not use the equipment without a locating cylinder.

Keep a distance of at least 5mm between the locator cylinder and the surface of the patient, so that it does not come into contact with non-biocompatible materials equipment.

Special handling precautions for the movable column version

During the handling of your X-ray equipment, be careful with moving parts indicated by the figure below, to avoid catching your fingers.

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**THIS EQUIPMENT CAN NOT BE
TILTED MORE THAN 5°.
RISK OF TUMBLING.**

Points to be observed prior to the first load application

Make sure that all the items concerning operation instructions have been completed.

- Operation (Seletronic)
- Operation (Eletronic)

Method used for measuring the X-ray tube voltage

Measurement is made directly in the high-voltage circuit by connecting a peak reading voltmeter between the tube anode and cathode by means of a resistive divider.

Method used for measuring the X-ray tube current

Measurement is made in the rectified high-voltage circuit. A continuous-current amperemeter is placed in series with the high-voltage transformer's secondary and ampoule, where the reading of the mean-value current is made.

Procedure in case equipment functioning is altered

In case the equipment presents any faults, firstly check whether they are listed under "Troubleshooting". If the problem cannot be solved, turn off the equipment, unplug it from the power source outlet and call for Dabi Atlante Authorized Technical Assistance.

Sensitivity to predictable environmental conditions in normal working situations

This equipment is not sensitive to magnetic, electric and electrostatic or pressure interferences as long as the cleaning, maintenance, transportation and operation recommendations found in this manual are observed.

Procedures in case of equipment disposal

When the equipment is put out of use, it must be disposed of at an appropriate site (according to local regulations), thus preventing its inadequate use and environmental contamination.

Transportation and storage conditions

The equipment must be transported and stored:

- Carefully in order to prevent falls and impact.
- By keeping the side indicated by the arrow facing upwards.
- By observing maximum piling limits as indicated on the packaging.
- Protected against moisture, rain, water sprinkling and wet ground.

Transportation /storage environmental conditions	
Transportation and storage environmental temperature range	0°C to + 55°C
Transportation and storage relative humidity range	0 % to 90 % (non-condensing)
Atmospheric pressure range	500hPa to 1060 hPa (375 mmHg to 795 mmHg)

After installation

Keep the equipment protected from sunlight and rain.

Environmental operation conditions	
Environmental temperature functioning range	+10°C to +35°C
Environmental temperature range recommended by Dabi Atlante	+21°C to + 26°C
Functioning relative humidity range	30 % to 75 % (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Environmental accommodation operation conditions (between operations)	
Environmental temperature accommodation range	+5°C to +45°C
Environmental temperature range recommended by Dabi Atlante	+15°C to +30°C
Functioning relative humidity range	30 % to 75 % (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Preparation prior to delivery

This equipment has been revised and approved according to the Dabi Atlante Quality Inspection and Control Program for better performance.

Dabi Atlante authorized service network

Installation and all maintenance service for Dabi Atlante equipment must be given by Dabi Atlante Authorized Technical Assistance, otherwise warranty will be compromised. Use the attached insert or access the website www.dabi.com.br to find the nearest Dabi Atlante Authorized Technical Assistant.

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**SPECIAL CARE FOR TRANSPORTATION
SEE MANUAL**

E

Fabricante
CNPJ: 55.979.736.0001/45
LIMITADOR DE FEIXE
Modelo: SPECTRO 70X

N° Série

Symbols

Use the icons below to identify your equipment's symbols.

	"Focal Point" - indicates the exact position of the radiation emission center.
	"Radiation" - indicates that the equipment emits ionizing radiation.
	"Shooting" - indicates the device for x-ray emission.
	This warns that, during transportation, the arms and tube head must be immobilized by a strap in the fully closed position, thus preventing the mobile unit from tilting.
	Indicates the tube head angle in relation to the horizon.
	Indicates the grounding terminals for protection.
	Caution - warning for consultation of enclosed documents.
	Type B - Applied Part
	"Alternate current" - indicates the equipment' input current type.
	"Continuous current" - indicates the equipment's output current type.
	"Fragile" Located on the packaging side. It determines that equipment must be carefully transported in order to prevent falls or shocks
	"Protect against moisture" Located on the packaging side. It determines that the equipment must be protected against any type of moisture during transportation and storage.
	"This side up" Located on the packaging side. It determines that the packaging must be handled by keeping the side shown by the arrow facing upwards.
	"Maximum piling" Located on the packaging side. It determines the maximum number of boxes that can be piled up during transportation and storage.
	"Temperature limit" It determines the temperature limits at which the equipment must be stored or transported.

Final recommendations

In order to ensure greater durability for your equipment, use only original Dabi Atlante replacement parts, which are in accordance with the technical standards and specifications required by the manufacturer. Also, use only the services provided by the Dabi Atlante Authorized Technical Assistance Network as it counts on trained assistants and specific tools for the appropriate maintenance of your equipment.

Equipment warranty

The "Warranty Terms" are found on the back of the certificate accompanying the product. Any alterations or concessions involving the warranty as well as pertaining conditions and authorizations, either oral or written, are definitely not permitted if not previously documented and approved of by Dabi Atlante.

Technical features

Models	Spectro 70X Seletronic Spectro 70X Eletronic Note: the standard model is supplied with a "movable column"
Optional	Wall mounted (without a column/base)

Classification according to Directive 93/42 Class II b

Classification according to standard NBR IEC 60601-1

Protection against electrical shock	TYPE-B APPLIED PART - Class I	
Protection against harmful penetration	Common level – equipment sealed without protection against water penetration	
Operation mode	Continuous operation in imminent action status for application of specific loads (intermittent load)	
Generator	Immerse in naftenic mineral isolating oil	
Protection level	Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide	
Cooling conditions	Automatically controlled cooling time. During the cooling time (15 seconds), the equipment remains locked for selection and activation	
Weights	60.9 Kg – movable column	Head: 8.4 kg Base: 37 kg
	30.9 kg – Wall-mounted	Timer case: 15.5 Kg Head: 8.4 kg Support/arm: 22.5 Kg

The x-ray emitting set is not pressurized.

Input Data

Mains Voltage	127/220V
Current	7.5 / 4.3 A~
Frequency	50 / 60 Hz
Power	1.20 kVA $\pm$ 10%
Number of phases	Monophase
Current type	AC (alternate current)
Voltage range	105 to 135V~ at 127V~ / 195 to 250V~ at 220V~
Fluctuation allowed for the voltage range	$\pm$ 3V~
Apparent resistance of power supply system	0,17 Ω for 127V~ 0,50 Ω for 220V~
Fuse	15A for 127V~ (normal fusion – glass 6.35 x 31.80 mm) 10A for 220V~ (normal fusion - glass – 6.35 x 31.80 mm)

Data on the High-Voltage Generator

X-ray tube nominal voltage	70 kVp $\pm$ 10%
Current	8 mA $\pm$ 20%
Power	0.41 kW (nominal)
X-ray tube model	D-082B-Toshiba or KL1-0.8-70-Kailong
Target material	Tungsten
Target angle reference axis	20° (D-082B-Toshiba) or 19° (KL1-0.8-70-Kailong) (in relation to anode)
Focal point	0.8x0.8 mm (positioned in relation to the reference axis) from the x-ray tube, according to IEC 60336

Overload Protectors		
127V~ voltage	15 A~ circuit breaker	
220V~ voltage	10 A~ circuit breaker	
Quality equivalence filtration	Toshiba Tube	Kailong Tube
Tube glass	1,41 mm Al	1,47 mm Al
Isolating oil (oil column 10.5mm)	0,67 mm Al	0,67 mm Al
Acrylic	0,25 mm Al	0,25 mm Al
Additional filter (aluminum)	0,70 mm Al	0,58 mm Al
Total filtration	3,03 mm Al	2,97 mm Al
Filtration with quality equivalence in 70 kV and 2.5 mm of aluminum with a semi-reducing layer – according to standard IEC 60522		
Data for application of a 0.3-second load		
Voltage	70 kVp ± 5kVp	
Current	8,0 mA ± 1,6 mA	
Power	0,54 kW	
Current time product	2.1 mAs.	
 - Continuous Thermal Maximum Dissipation (CTMD) = 2 mA x 70 kV = 140W - The field specified for radiation time conformity ranges from 0.3 to 3,2 seconds.		

Load application parameter values related to leakage radiation:

According to standard NBR IEC 60601-1-3, the values for leakage radiation do not exceed the limit of 0.25 mGy in one hour to 70 kV parameters.

Conformity Statement

Dabi Atlante guarantees that the Spectro 70X x-ray machine, as concerns all models cited in this manual, conform to standard NBR IEC 60601-1:1997

Radiation source set X SPECTRO 70X, NBR IEC 60601-2-28:2001

Troubleshooting

Faults	Causes	Solutions
Keys will not work	Equipment in unplugged	Plug it in
	Fuse is burned	Replace fuse
	Master switch is off	Turn switch on
	Circuit breaker is off	Turn on circuit breaker
	No power supply	Wait until supply is normalized
A semicircle appears on radiographs	Locating cylinder incorrectly positioned	Apply the parallelism radiographic technique by using the auxiliary lines of the collimating cylinder
Radiographs show a dark strip	Dark room	Light has penetrated the dark room. Prevent sunlight or excessive illumination from entering the dark room
Control shows a dark strip	Faulty electronic circuit	Unplug the equipment and call for technical assistance
Completely dark radiographs action	Excessive exposure time	Adjust exposure time (refer to "Radiographic Techniques" in this manual)
	Excessive developing time	In hot weather, the developer has quicker
	Developer Note: Kodak developer does not use mixtures	Maybe the developer may contain insufficient mixture (water/ developer ratio – 3:1)

Troubleshooting

Faults	Causes	Solutions
Keys will not work	Equipment in unlugged	Plug it in
	Master switch is off	Turn switch on
	Circuit breaker is off	Turn on circuit breaker
	No power supply	Wait until supply is normalized
Over or underexposed film	The kV adjustment value is not adequate for the patient's size. Possible error during film development	Check the development system or repeat the exam varying the KV value. If this procedure fails, call for technical assistance
Darkened radiograph	Definition failure on the whole film	Check if the film's best-before date has expired and if storage conditions are being followed. Check if the light bulb in the dark room is appropriate and if light is penetrating the dark room
Radiograph with strips	The film may have been consecutively rubbed in the packaging case and accumulated electrostatic charges, or the reinforcement screens were not treated by the anti-static solution	Avoid rubbing film in the case and clean the screens with anti-static solution
Radiograph with blurred details and good contrast	Film is not correctly pressed between the screens	Replace cassette
Radiograph is completely blank	No x-ray emission or development process is irregular	Check x-ray passage by means of a sound and light signal. In case error persists, check the development process. If problem is still not solved, call for technical assistance

Appendices

The Dental X-Ray Spectro 70X models Seletronic and Eletronic must be used in an electromagnetic environment, such as that specified below. The client or user of the Dental X-Ray Spectro 70X models Seletronic and Eletronic must ensure that it is used in such type of environment.

Directives and manufacturer's statement – Electromagnetic emissions

Assays and emissions	Compliance	Electromagnetic environment - directives
RF emissions ABNT NBR IEC CISPR 11	Group 1	The Dental X-Ray Spectro 70X models Seletronic and Eletronic uses RF power only for its internal functions. However, its RF emissions are very low and it is not likely that they will cause any interference with nearby electronic equipment.
RF emissions ABNT NBR IEC CISPR 11	Class A	The Dental X-Ray Spectro 70X models Seletronic and Eletronic is suitable for use in all types of facilities, including residential facilities and those directly connected to the public service network of low-voltage power supply providing service to residential buildings
Harmonics emissions IEC 61000-3-2	Class A	
Emissions due to voltage fluctuation / scintillation IEC 61000-3-3	Compliant	
RF emissions CISPR 15	Compliant	The Dental X-Ray Spectro 70X models Seletronic and Eletronic is not appropriate for interconnection with another piece of equipment.

Directives and manufacturer's statement – Electromagnetic immunity

The Dental X-Ray Spectro 70X models Seletronic and Eletronic must be used in an electromagnetic environment, such as that specified below. The client or user of the Dental X-Ray Spectro 70X models Seletronic and Eletronic must ensure that it is used in such type of environment.

Immunity assays	Assay level of ABNT NBR IEC 60601	Compliance level	Electromagnetic environment - directives
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV through contact ± 8 kV through air	± 6 kV through contact ± 8 kV through air	Floors must be made of wood, concrete or ceramic. If they are finished with synthetic material, relative moisture should be of at least 30%.
Fast electric transitories/pulse train (Burst) IEC 61000-4-4	± 2 kV on the mains lines ± 1 kV on the input/output lines	± 2 kV on the mains lines ± 1 kV on the input/output lines	Power supply quality must be that of a typical hospital or commercial environment
Outbreaks IEC 61000 - 4 - 5	± 1 kV differential mode ± 2 kV regular mode	± 1 kV differential mode ± 2 kV regular mode	Power supply quality must be that of a typical hospital or commercial environment
Voltage failure, short interruptions and voltage variations in the mains input line IEC 61000-4-11	<5% U_T (>95% of voltage failure in U_T) for 0.5 cycle <40% U_T (>60% of voltage failure in U_T) for 5 cycles <70% U_T (>30% of voltage failure in U_T) for 2.5 cycles <95% U_T (>5% of voltage failure in U_T) for 5 cycles	<5% U_T (>95% of voltage failure in U_T) for 0.5 cycle <40% U_T (>60% of voltage failure in U_T) for 5 cycles <70% U_T (>30% of voltage failure in U_T) for 2.5 cycles <95% U_T (>5% of voltage failure in U_T) for 5 cycles	Power supply quality must be that of a typical hospital or commercial environment. If the user of the Dental X-Ray Spectro 70X models Seletronic and Eletronic requires continuing operation during a power failure, it is recommended that the Dental X-Ray Spectro 70X models Seletronic and Eletronic be supplied by an uninterrupted power source or a battery.
Magnetic field in the supply frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields with frequency in the supply must be in characteristic levels of a typical site in a hospital or commercial environment
			Portable and mobile RF communication equipment must not be used near the Dental X-Ray Spectro 70X models Seletronic and Eletronic (including cables) with the recommended separation distance calculated from the equation applicable to the transmitter's frequency. Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz up to 2,5 GHz	3 Vrms	$d=1,2\sqrt{P}$

Radiated RF IEC 61000-4-3	3 Vrms 80 kHz up to 2,5 GHz	3 V/ms	$d=1,2\sqrt{P}$ - 80 MHz up to 800 MHz $d=1,2\sqrt{P}$ - 800 MHz up to 2,5 GHz Where P is the transmitter's maximum nominal output power in watts (W), according to the transmitter's manufacturer, and d is the recommended separation distance in meters (m). It is recommended that the field intensity established by the RF transmitter, as determined by an electromagnetic inspection on the site ^a , should be lower than the compliance level in each frequency range ^b . Interference may occur around equipment marked by the following symbol: 
------------------------------	--------------------------------	--------	---

Note 1: U_T is the AC supply voltage prior to the application of the assay level.

Note 2: In 80 MHz 800 MHz, the highest frequency range is applied.

Note 3: These directives may not be applicable in all situations. Electromagnetic propagation is affected by the adsorption and reflection of structures, objects and people.

^a The field intensities established by fixed transmitters such as base radio stations, telephone (cellular/cordless), terrestrial mobile radios, amateur radio, AM and FM radio transmission and TV transmission cannot be accurately predicted. In order to evaluate the electromagnetic environment due to fixed RF transmitters, it is recommended to consider an electromagnetic inspection of the site. If the field intensity measurement on the site where the Dental X-Ray Spectro 70X models Seletronic and Eletronic is used exceeds the RF compliance level applicable above, it is recommended to check the Dental X-Ray Spectro 70X models Seletronic and Eletronic in order to observe if operation is normal. In case abnormal performance is observed, additional procedures may be necessary, such as reorientation or replacement of the Dental X-Ray Spectro 70X models Seletronic and Eletronic.

^b It is recommended that field intensity should be smaller than 3V/m over the frequency range of 150 kHz to MHz.

Recommended separation distance between portable and mobile RF communication equipment and the Dental X-Ray Spectro 70X models Seletronic and Eletronic

The X-Ray Spectro 70X models Seletronic and Eletronic has been designed for use in electromagnetic environments where radiated RF disturbances are controlled. The client or user of the X-Ray Spectro 70X models Seletronic and Eletronic can help prevent electromagnetic interference by keeping a minimum distance between portable and mobile RF communication equipment (transmitters) and the X-Ray Spectro 70X models Seletronic and Eletronic as recommended below, according to the maximum output power of the communication equipment.

Transmitter maximum nominal output power (W)	Separation distance according to transmitter's frequency (m)		
	150 kHz to 80 MHz $d=1,2\sqrt{P}$	80 MHz to 800 MHz $d=1,2\sqrt{P}$	800 MHz to 2,5 GHz $d=2,3\sqrt{P}$
0,001	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters whose maximum nominal output power is not listed above, the recommended separation distance d in meters (m) can be determined by using the equation applicable for the transmitter frequency, where P is the transmitter's maximum nominal output power in watts (W), according to information provided by the transmitter's manufacturer.

Note 1: In 80 MHz 800 MHz, the highest frequency range is applied.

Note 2: These directives may not be applicable in all situations. Electromagnetic propagation is affected by the adsorption and reflection of structures, objects and people

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MPC International S.A - 26, Boulevard Royal - L-2449 - Luxembourg